

SCLERAL BUCKLING PRODUCTS


Packaging: **Box/5**
Sterile: **Yes**

Our extensive range of scleral buckling products consists of both solid and sponge implants and includes circling bands, tires, silicone strips, sleeves and silicone sponges.


SOLID SILICONE IMPLANTS

Silicone Strips

Silicone Strips can be used alone, with wide grooved tires, or when 'trap door' procedures are performed.




Article No: 92-09
Style 41, 3.5mm Silicone Strip
41




Article No: 92-10
Style 42, 4.0mm Silicone Strip
42

Round Silicone Sleeves

Sleeves are used to secure circling bands.




Article No: 92-13
Style 70, Silicone Sleeve
70




Article No: 92-14
Style 270, Silicone Sleeve
270

Silicone Circling Bands

To hold the accessory implant in place and to maintain optimal pressure to form a permanent buckle without scleral erosion.




Article No: 92-01
Style 40, 2.0mm Circling Band
40




Article No: 92-02
Style 240, 2.5mm Circling Band
240

Silicone Tires


Tire-shaped implants are designed for use under circling bands for wide scleral beds, for breaks near the ora serrata, multiple breaks, and for high scleral buckles.




Article No: 92-15
Style 277, Concave Tire
277




Article No: 92-16
Style 279, Concave Tire
279



Article No: 92-17
Style 287, Convex Tire
287



Article No: 92-18
Style 276, Asymmetrical Tire
276




Article No: 92-22
Style 286, Convex Tire
286

SCLERAL BUCKLING PRODUCTS


SILICONE SPONGES

Silicone Sponges are offered in many shapes that can be used for segmental buckling procedures.


Round Silicone Sponges



Article No: 92-50-3R
Style 503, 3.0mm Round
Sponge - 503




Article No: 92-50-4R
Style 504, 4.0mm Round
Sponge - 504




Article No: 92-50-5R
Style 505, 5.0mm Round
Sponge - 505

Oval Silicone Sponges



Article No: 92-50-5
Style 506, 3.0mm x 5.0mm
Oval Sponge - 506




Article No: 92-50-7.5
Style 507, 5.5mm x 7.5mm
Oval Sponge - 507

ACCESSORIES

Meridional Implants

Gives additional buckling in the meridional direction.



Article No: 92-27
Style 112, Meridional Implant
12.0mm -112

For additional information please contact:
D.O.R.C. Dutch Ophthalmic Research Center (International) B.V.
Scheijdelveweg 2, 3214 VN Zuidland, The Netherlands

Phone: +31 181 45 80 80 | Fax: +31 181 45 80 90 | E-mail: sales@dorc.eu